# Serie 1000

**RU|MED**®  
Rubarth Apparate GmbH

## Licht-Thermostate Light Thermostats

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## RUMED Light Thermostats - Limitless Individuality

RUMED Light Thermostats of Series 1000 - a variable programme with a multitude of options for special requests and for limitless individuality. The standard series can be adapted flexibly to the specific requirements of the user by adding a special lighting equipment and the corresponding accessories: programme control and humidity.

RUMED Light Thermostats are an efficient alternative to the Cooled Incubator, if high lighting intensities for temperature application are required.

Advantage: flexibility and highest precision.

### Various Application Ranges

For tasks in research and quality control of industry and science in the ranges:

- Plant research
- in vitro propagation
- Seed Tests and Seed Research
- Zoology and Botany
- Photo stability test

### General view of the most important features

- Test room volume 210 l to 1.340 l
- Test room temperature -20°C to +50°C
- High operating convenience
- Intuitive operation
- Solid and longevous
- Easy to repair and to maintain
- Possibility of calibration and validation
- Upgrade capability

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<b>The Appliance - Big Volume on a Small Floor Space</b>	The high test room is a main feature for favourable proportion of all RUMED precision appliances. Advantage: big volume and small floor space. The test room door can be widely opened to reach the whole cross-section of the interior space. A PVC magnetic seal surrounding the door frame seals the door reliably. For appliances with test room temperatures	below 0 °C a controlled door frame heating avoids formation of condensation water in the area of the door seal. The exterior housing of all models is fabricated in stainless steel execution. A particularly efficient insulation is achieved by high-quality CFC-free polyurethane foam, so that the preset temperature can be kept constant with minimum expenditure of energy.
<b>The Test Room - Easy-Care and Variable</b>	The test room of all models consists of stainless steel. The entire test room, as well as the shelves, are easy to clean due to their smooth surface and their rounded edges and corners. The shelves consist of plastic-laminated wire. They can be delivered in stainless steel execution upon request. 4 different test room variants are available: Appliances with lighting from above in	one level for high growth height, appliances with lighting from above in several levels for low growth height, appliances with lighting from both sides and possibility of charging in several levels, and appliances with special air ducts for condensate-free in vitro propagation.
<b>The Heating - Continuously Adjustable and Resistant to Wear</b>	Coasting of temperature is almost avoided due to the low mass of the electric resistance heater. The heater is arranged directly in front of the recirculating air blower. The continuous and contactless control by means of a solid state relay allows an extremely direct and precise control.	Advantage: High temperature accuracy and minimum wear.

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**The Cooling -  
Energy-Saving and CFC-Free**


The air-cooled refrigerating machine is working with CFC-free refrigerant R134a (in case of temperature ranges below 0°C R404a). The compressor is running quietly and largely free from vibrations. The efficiency of the heat-exchanger is optimal. The total cooling circuit is working with solenoid valve-bypass technique ensuring that the compressor will only be disconnected, if cooling capacity has not been required for a prolonged period.

Result: Reliability, high operating safety and durability.

**The Temperature -  
Precise and Constant**

A precision appliance is expected to be reliable and precise. Our appliances guarantee high temperature accuracy in time and in space according to the

DIN standard 58945 due to the continuously running fan, which can also be delivered with speed control for most of the units .

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**1001 with Options**  
 Temperature range -20 to +50°C  
 Test room volume 210 litres  
 Lighting from above

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**1101 with Options**  
 Temperature range -20 to +50°C  
 Test room volume 350 litres  
 Lighting from above

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**1201 with Options**  
 Temperature range -20 to +50°C  
 Test room volume 350 litres  
 Lighting from both sides

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**1301 with Options**  
 Temperature range -20 to +50°C  
 Test room volume 530 litres  
 Lighting from both sides

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
**Dialogue with the User -**

Plain text menus with selection by luminous bars ensure easy operation. Two languages are standard. Languages like German, English, French, Italian and Polish are at choice.

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
**Structured Programming -**

Tables are facilitating the programming. Numerous editing functions are offered.

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

**Easy Verification -**

Plausibility check with indication of minimum and maximum setpoint.

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**Constant observation -**  
of the programme execution by means of graphic display; display of the actual moment; possibility to scroll through the different profiles and days.

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**Calibration and adjustment -**  
of the temperature and humidity sensors by integrated supporting point correction - The determined correction values for the limit ranges and the values at the three freely selectable supporting points can be entered easily into a clearly arranged table - They can be read at any time.

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
**Appliance logbook -**  
with 99 possible entries - priority controlled FIFO function - Selective acknowledgement or deletion of messages and alarms, which can be read-out via interface.

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
Individual units can be supervised by the standard interface RS232. A serial computer interface RS485 is available as option. A RS232/RS485 interface converter allows the management of up to 32 units each.

The optional software package allows:

- Data logging (monitoring) of the data being relevant to the process
- Programme management (reading, writing, transmission, storage etc.)
- Log of messages and alarms with time stamp (logbook function)
- Display and print-out of diagram evaluations of the recorded data

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
**Standard Unit  
Single Setpoint Control**


**Standardgerät  
Photoperiode**

**Standard Unit  
Photoperiod**


**Option Programmregelung  
70-Stufen Programm**

**Option Programme Control  
70-Stage Programme**

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---

## The Lighting - The Right Variant for Every Application

For plants with a large growth height, units with lighting from above are available. In this case, only one lighting level is recommended. For plants requiring a particularly high luminous intensity, the unit can be equipped with metal halide vapour discharge lamps. Depending on accessories and model, maximum luminance values of up to  $800 \mu\text{E}/\text{m}^2$  (equivalent to approx. 40.000 Lux) can be attained.

Appliances with lighting from above in several levels are the best choice for plants with low growth height. Optimum economy of space and economy of energy are the main features of these appliances.

The universal cabinets with lighting from both sides are suitable for small specimens. These appliances offer a free and flexible arrangement of the test room. Additional shelves ensure a maximum of light on test room surface and a minimum floor space being required for the appliance at the same time.

The fourth appliance execution has been particularly developed for "in vitro" cultivation. The appliance has been designed with a special air conduction and lighting from above on several levels.

All these appliances are prepared for selection of the luminous intensity in 2 or more stages. If the standard appliances are equipped with the Control2000, a simple photoperiod can be run.

In combination with the option programme control and the option continuous control of luminous intensity (only in case of lighting with fluorescent lamps), complex temperature and lighting profiles can be run.

In case of appliances with lighting equipment being mounted outside the test room, a double glazing effects an efficient thermal insulation between the lighting equipment and the test room. This is energy saving, since the test room temperature can be kept constant with low energy consumption. At the same time, high-efficiency reflectors ensure an optimum luminous efficiency.

Lighting is available in different light colours, the spectra of which are shown in the graphs on page 11. Apart from the shown spectra, further light colours are available for appliances with fluorescent lamps. (Please contact us!)

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

**Option:  
Humidification by Evaporation**

Two versions are available for application with humidification: a low-cost version for simple tasks and a version for tasks requiring highest precision. For the low-cost version, an evaporation humidifier is installed in the test room. Control of the relative humidity is realized by the PID-controller being integrated in the CONTROL2000. This execution is always recommendable, if high air humidity values are to be kept (higher than the ambient air humidity), and if high precision is not requested.

The sensor is a temperature-compensated, capacitive sensor, which can be calibrated.

The water supply is either realized automatically by direct connection to the drinking water supply or by periodical manual refilling of drinking water.

**Option:  
Humidification by Ultrasonic  
Atomization**

The version for application with high precision requirements is working with ultrasonic humidification and connection to demineralized water. The humidity control is also effected by the electronic PID-controller being integrated in the Control2000. This method of humidification combined with the option dehumidification allows preselection of the desired humidity according to the diagram in the margin. The option ultrasonic humidification effects the optional dehumidification by condensation at a specially developed plate heat-exchanger of stainless steel.

The response characteristic of this method is substantially improved (quick change of humidity) and substantially lower air humidity values are achieved. For high requirements, only the combination of humidification and dehumidification is recommendable, in order to get a real independence of ambient conditions.




The sensor is a temperature-compensated, capacitive sensor, which can be calibrated. The humidity sensor is equipped with a bayonet catch. Thus, it can be replaced without any tool.

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## Spectral Distribution of Radiation of the Lighting Equipment

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


**1501 with Options**  
 Temperature range 0 to +50°C  
 Test room volume 800 litres  
 3 Levels with lighting from above

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---

**For Appliance and Specimen -  
Double Safety**

Superheating and supercooling is no more existing, since failures are detected by additional, independent excess temperature and insufficient temperature sensors, and the Control2000 disconnects the corresponding units reliably. In addition to the independent appliance fuse TS (DIN 12880, class 1) existing in all appliances, the appliance can be equipped

with an optional, independent temperature selector according to DIN 12880, class 3.3. This represents an inquiry of safety of the integrated temperature control of the Control2000.

**Option:  
Recirculating air as requested -  
The right wind force for everyone**

The option continuously adjustable recirculating air combined with the programme control allows different intensity values of recirculating air. The linearization of the controllable range can be adapted to the own measuring values by means of a supporting point correction.

Advantage:  
Flexible adaptation to individual requests.

**Option:  
Fault Message Contact -  
Inquiry of Safety**

The appliance can be equipped with a potential-free fault message contact. Apart from local monitoring, the appliance can be connected to a centralized control room by means of that contact. Except from power failure, configuration in several steps can be effected: if messages shall only be displayed locally, and which alarm messages shall be displayed in the external control room.

Advantage:  
Advanced safety for unit and specimen due to early recognition and premature possibility to react.

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**In Vitro Propagation -  
Free From Condensate is Healthy**

The special air conduction through air ducts avoids condensation in the culture vessels during in vitro cultivation. A continuous airflow from below to above is guided around the cultures being placed on a grating above the air ducts. Since formation of condensate drops is avoided, the cultures will not be contaminated by dripping water. Further, the burning glass effect of the lighting is avoided and the free view on the cultures for examination is always ensured.

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**Rotational Stand -  
Motion is Required**

A powered rotational stand ensures uniform distribution of the light on the algae cultures. The individual levels consist of acrylic glass to minimize shading. The 6 levels are equipped with 8 supporting devices, each, for the culture vessels. The rotational stand can be equipped with an adjustable continuous drive, as well as with an interval drive, which is avoiding the sedimentation. An optional gas supply device ensures a separate supply of the individual cultures. Admission is effected centralized through the shaft of the rotational stand. The desired quantity of gas for the individual cultures can be supplied through 48 needle valves being individually adjustable.

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The figured units are standard units of the stated types with options. The indicated temperature range is based on an ambient temperature of +23°C without heat emission of the specimen. All dimensions, connected load rates and indications of weight are approximated values and are valid for the standard version. Technical and optical changes due technical progress are subject to change without notice.

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Optional Equipment		
Movable design by 4 swivelling guide rolls.		1900
Switchable socket for the operation of electric units in the test room (f. ex. shaker). Protected by a 4A safety cut-out.		1910
Tubular port for measuring lines into the test room (20 mm diameter).		1915
Tubular port for measuring lines into the test room (45 mm diameter).		1916
Speed control of the fan. Reduces the air velocity for sensitive specimens, avoids drying. Not available for model 1601, 1701 and 1801!		1925
Separate Working Hour Meter, which can be reset by key switch.		1930
Temperature Selector Switch TWW being independent of the controller, according to DIN 12880, class 3.3, as further protection against controller malfunction.		1950
Ion exchanger for direct connection to the drinking water supply, if demineralized tap water is not available. Only for units with the option ultrasonic humidification.		1970
2 Spare cartridges for ion exchanger. Worn cartridges are taken back free of charge. The cartridges will be recycled.		1971
Canister for supply of the ultrasonic humidifier for locations without water supply.		1975
Pump for Water Supply from Supply Vessel.		1978
Heated condensation water bowl for automatic evaporation of the accumulated condensation water.		1979
Programme controller for daily-cycles, weekly-cycles, realtime and processtime programmes. Programme with 70 independent steps in every profile (f. ex. temperature, humidity, light, socket etc.)		1981
Programme controller same as above but 4 programme memory locations. Many additional possibilities (cycles of 2 or even more programmes).		1982
Serial digital computer interface RS 485 for integration of the unit into a software monitoring system.		1986
Combined Voltage/Current Output (0-10 V, 0-20 mA, 4-20 mA) for recording of the test room temperature by means of an external recorder.		1990
Combined Voltage/Current Output (0-10 V, 0-20 mA, 4-20 mA) for recording of the test room humidity by means of an external recorder.		1991
1-Channel recorder for documentation of the temperature course according to the GLP and GMP regulations. Will be mounted in the switchboard.		1992
2-Channel recorder for documentation of the temperature and the humidity course according to the GLP and GMP regulations. Will be mounted in the switchboard.		1993
Potential-free contact, which can be freely configured, for failure alarm in a control room		1995
Software for centralized monitoring, storing of recipes and data-logging.		6000
Calibration of one unit.		6700
Validation of one unit.		6800

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Test and Simulation Appliances  
for Research, Quality Control  
and Production

DE-S-133 00001