Serie 5000

RU MED®
Rubarth Apparate GmbH

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RUMED Series 5000 - Systematic Seed Testing

Rumed seed testing equipment of series 5000 - an extensive programme for systematic seed testing according to the ISTA standard - from sample division to determination of germination.

The high-quality execution, as well as the simple handling and the thought-out construction are convincing.

Various Application Ranges
For tasks in research and quality control of seed testing and seed production

General View of the Most Important Features:

Rumed germinators are the professional alternative to cooled incubators, wherever high capability of reproduction of the results with low lighting intensities are required.

Advantage: Reliability and durability, as well as highest precision at the same time.

- **Sample division**
- **Purity examination**
- **Determination of germination**

- **Effective surface for 120 to 180 samples**
- **Water temperature from +5 °C to +40 °C**
- **High operating convenience, ease of operation**
- **Easy cleaning**
- **Easy to repair and to maintain**
- **Recycling of old units**

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The Jacobsen Method

The Jacobsen germinator mainly consists of a germination plate being temperature-conditioned by means of the water basin below. The water bath is equipped with an automatic temperature control.

The germination spirals being equipped with a paper wick and a paper substrate are placed on the germination plate. The wick is being led through slots in the germination plate and reaches into the water bath below, thus supplying the required humidity and the desired temperature to the paper substrate.

The two circular filter papers are covered with a transparent cover dome to provide the air humidity being required for the germination. A little hole in the upper end of that dome ensures sufficient supply of fresh air and minimum evaporation at the same time.

Units being executed with active cooling allow day-night temperature alternation, as well as any temperature profile.

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The Rodewald Method

Rodewald germinators are mainly consisting of a tub being filled with silica sand, which is hanging in a temperature conditioned water basin for temperature conditioning and moistening.

The Rodewald method is particularly suitable for vegetable or forest seed due to the adjustable low humidity and the inferior tendency to fungal growth. The seed is deposited on filter papers on the sand.

A water channel being integrated in the sand insertion tub with adjustable level regulation gives moisture to the seed via wicks in the sand and the capillary effect of the sterilized crystal silica sand.

An acrylic glass cover avoids drying-out of the sand. The glass cover can be steplessly opened for ventilation. A surrounding groove avoids dripping of the condensate water into the sand.

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5401
Rodewald Germination Table
with opened acrylic glass cover

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

Dialogue with the User -
Plain text menus with selection by luminous
bars ensure easy operation. Two languages
are standard. Languages like German,
English, French, Italian and Polish are at
choice.

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
Structured Programming -
Tables are facilitating the programming.
Numerous editing functions are offered.

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Easy Verification -
Plausibility check with indication of minimum
and maximum setpoint.

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**The Water Tub -
Guarantees good temperature-
conditioning**

The one-piece water tub consists of stainless steel sheet, material No. AISI 316Ti. The thermal insulation of the exterior tub surface reduces the exterior temperature influence to a minimum.

Balancing of evaporation loss is effected by means of an opto-electronic level regulator keeping the water level automatically on a constant level. A water overflow guarantees highest possible safety.


The water supply installation of the unit is executed according to the German regulations for gas and water installation DVGW. The unit can be connected permanently to the drinking water supply. A connection to the sewage or a floor drain is required.

Changing of the water is effected by actuation of a drain valve. Filling is effected unattended due to the automatic filling time control.

**The Connections -
A Good Connection is the Main Thing**

All electric and sensory connections are waterproof plug-in screwed connections being easy to maintain.

All water connections are executed as quick-fitting pipe unions. Cleaning of the pump can be effected without any problems. It can be removed and installed without any tools in only 1 minute.

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Temperature Plus Lighting

Even the standard single setpoint control in germinators with lighting allow running of a photo period due the real-time clock being integrated in all units. Combined temperature-lighting profiles require the execution with "Programme Control".

All Jacobsen germinators are prepared for unproblematic retrofit of lighting.

Specimen Sensor

All units are equipped with a specimen sensor. The temperature measured by the specimen sensor is displayed at the


Display as well as the temperature of the water bath.

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
5101
Jacobsen Germination Table
with lighting

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
Individual units can be supervised by the standard interface RS232. A serial computer interface RS485 is available as option. A RS232/RS485 interface converter allows the management of up to 32 units.

A software package for supervision and for central data logging is also available.

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Soil Divider for Samples

The soil divider for samples being completely executed in stainless steel makes the reliable and representative sampling of partial samples dead easy. The seed is just to be placed into the opened cover and it is to be distributed more or less uniformly. Due to the special construction with 18 channels, the seed will be separated into two partial samples of the same size as soon as the cover is closed.

After the separation, the partial samples can be withdrawn easily from the two lateral drawers.

The void-free construction ensures easy cleaning without any residues. Due to the electropolished surfaces sticking of residues is practically made impossible.

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Vacuum Seed Counter

This extremely practical device facilitates the tiring and time-consuming counting and uniform depositing of seed for germination. It is particularly suitable for seed being regularly shaped and relatively even, such as grain, brassica and trifolium sorts. The seed counter consists of three main components: a vacuum system with connection hose, a number of counting heads corresponding to the different seed sorts and a valve to

dissolve the vacuum. The vacuum capacity can be adjusted by a potentiometer. The counting heads with 100, 50 or 25 bore holes are some smaller than the paper substrate and are equipped with an edge to prevent the seed from rolling off. The diameter of the holes is adapted to the size of the seed and the suction capacity of the vacuum.

Seed Blower

The seed blower serves for quick and easy separation of light and heavy seed. The seed is filled into a drawer of the acrylic glass cylinder. A fine-meshed sieve closes the upper end of the cylinder. A speed-controlled blower in the cylinder generates a constant upward air column. Depending on its size and weight, the seed will fall into one of the 3 collecting vessels being installed at the cylinder.

The light seed will be deposited in the upper collecting vessel, whereas the heavy seed will be deposited in the lower collecting vessel. A further sieve at the bottom of the cylinder avoids penetration of seed into the blower. After disconnection of the blower, the sieve contains also the deposit, consisting of the residues and eventually of small stones. The collecting vessels are equipped with a snap buckle clip allowing easy removal.

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
5851
Seed Blower
 For Separation of Light and
 Heavy Seed

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5801
Vacuum Seed Counter
 Facilitates Counting and
 Deposit of Seed

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5881
Soil Divider acc. to ISTA
 Cereal Sample, max. 2 kg

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- standard
- not available

The figured units are standard units of the stated types with options. The indicated temperature range is based on an ambient temperature of +23 °C. All dimensions, connected load rates and indications of weight are approximate values and are valid for the standard version. Technical and optical changes due to technical progress are subject to change without notice.

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Test and Simulation Appliances
for Research, Quality Control
and Production

DE-S-133 0000 1